

SAMPLE SUBMISSION GUIDE

Cloning Service

Sending samples according to the requirements below helps us to do our job better and provides you with accurate results.

Sample Data and Formats

- Please use our online forms for cloning service orders and quote requests
- We will calculate the price of the project based on your request and send you a Quotation. The Quotation summarizes your request and specifies the price of the service.

Sample Requirements

- Please send us all samples at a concentration of 200ng/ μ l
- The total sample volume should be minimum 10 μ l
- Please do not perform a restriction digestion of your vector beforehand
- The sample should be sent in doubled distilled water or TE buffer
- Please send your samples in 1.5ml tubes
- Labelling of the tubes:
 - Sample Name as indicated in the request form
 - Concentration of the sample
 - Volume of the sample

Shipping

- Please Mention the Order ID
- Please use a padded envelope when sending the vector to ensure that it arrives safely
- Ship the samples liquid at ambient temperature
- Shipping address:

Gene Synthesis - Lab
Eurofins Genomics
c/o Sarah Ebner
Anzinger Str. 7a
D-85560 Ebersberg

Contact

Please do not hesitate to contact us in case of any questions!

Email: genes-eu@eurofins.com; Phone: +49 7531 816 068